

COMPARATIVE STUDIES

Chelating Agarose Beads



RESINS


Experiment:

Unpurified extract containing Glutaryl acylase (6xHis) was tested under the same conditions with different **NICKEL** charged chelating beads. The SDS-PAGE shows the eluted fraction in all the resins.

Comparison:

- 1.- Low Molecular Weight markers (LMW)
- 2.- Glutaryl acylase (6x His) extract.
- 3.- [ABT](#) Cat. N°: 6BCL-QLNi
- 4.- *Competitor A*

	COMPANY	AVAILABLE	PURIFICATION OF <i>[Glu-6xHis]</i>
3	ABT	"READY to use" version	ONE SINGLE STEP PURIFICATION. HIGH YIELD OF PURIFIED PROTEIN
4	A	User has to bind Nickel before the purification procedure	Fused Protein very impure. Very Low Yield of Purified Protein




RESINS

Experiment:

Unpurified extract containing Fuculose 1-aldolase (6xHis) was tested under the same conditions with different **NICKEL** charged chelating beads. The SDS-PAGE shows the eluted fraction in all the resins.

Comparison:

- 1.- Low Molecular Weight markers (LMW)
- 2.- F1PA (6x His) extract
- 3.- *ABT* Cat. N°: 6BCL-QHNi
- 4.- *Competitor S*



	COMPANY	AVAILABLE	PURIFICATION OF <i>[F1PA-6xHis]</i>
3	ABT	"READY to use" version	ONE SINGLE STEP PURIFICATION. VERY HIGH YIELD OF PURIFIED PROTEIN
4	S	User has to bind Nickel to the resin	Fused Protein very impure. Very Low Yield of Purified Protein


RESINS

Experiment:

Unpurified extract containing Glutaryl acylase (6xHis) was tested under the same conditions to different COBALT charged chelating beads. The SDS-PAGE shows the eluted fraction in all the resins.

Comparison:

- 1.- Low Molecular Weight markers (LMW)
- 2- Glutaryl acilase (6x His) extract
- 3.- [ABT](#) Cat. N°: 6BCL-QHCo
- 4.- *Competitor S (Product 1)*
- 5- *Competitor S (Product 2)*



	COMPANY	AVAILABLE	PURIFICATION OF <i>[Glu-6xHis]</i>
3	ABT	"READY to use" version	HIGH YIELD & VERY GOOD SELECTIVITY OF PURIFIED PROTEIN
4	S	"READY to use" version	Very Low Yield of Purified Protein
5	S	User has to bind Cobalt to the resin	Very Low Yield of Purified Protein


RESINS

Experiment:

Unpurified extract containing Glutaryl acylase (6xHis) was tested under the same conditions with different COPPER charged chelating beads. The SDS-PAGE shows the eluted fraction in all the resins.

Comparison:

- 1.- Low Molecular Weight markers (LMW)
- 2- Glutaryl acylase (6x His) extract
- 3.- *ABT* Cat. N°: 6BCL-QLCu
- 4.- *Competitor A*



	COMPANY	AVAILABLE	PURIFICATION OF <i>[Glu-6xHis]</i>
3	ABT	"READY to use" version	VERY HIGH YIELD & VERY GOOD SELECTIVITY OF PURIFIED PROTEIN
4	A	"READY to use" version	Very Low Yield of and a significant amount of background proteins are also retained

In spite of the generally low selectivity of Copper, ABT's COPPER charged resin allows a HIGH DEGREE OF PURIFICATION OF THE TARGET-PROTEIN.


RESINS

Experiment:

Unpurified extract containing Glutaryl acylase (6xHis) was tested under the same conditions with different COPPER charged chelating beads. The SDS-PAGE shows the eluted fraction in all the resins.

Comparison:

- 1.- *Competitor A*
- 2.- *Competitor S*
- 3.- *ABT* Cat. N°: 6BCL-QLCu



	COMPANY	AVAILABLE	PURIFICATION OF <i>[Glu-6xHis]</i>
1	A	User has to bind Copper to the resin	Very Low Yield and a significant amount of background proteins are also retained
2	S	User has to bind Copper to the resin	Fused Protein impure
3	ABT	"READY to use" version	ONE SINGLE STEP PURIFICATION