



## TECHNICAL SPECIFICATIONS

### CHELATING Agarose Pre-packed XL Columns

BEAD (Geometry, size)	Spherical, Standard: ~ 50-150 µm					
CROSS-LINKED	Yes					
AUTOCLAVABLE	121°, 30'					
AGAROSE %	6% agarose					
ACTIVATING AGENT	Epichlorohydrin					
MATRIX	Stable in all commonly used reagents <sup>(1)</sup>					
BED VOLUME	5 ml					
COLUMN MATERIAL	Polypropylene column and polyethylene frit					
PRODUCT NAME	His-XL COLUMN HIGH Density			His-XL COLUMN LOW Density		
	NICKEL	ZINC	COBALT	NICKEL	ZINC	COPPER
BINDING / LOADING CAPACITY (µmol Me <sup>2+</sup> /ml gel)	20 - 40			5 - 20		
CAT. No.	6BCL-QHNI-C5	6BCL-QHZn-C5	6BCL-QHCø-C5	6BCL-QLNI-C5	6BCL-QLZn-C5	6BCL-QLCu-C5
ANTIMICROBIAL AGENT	20% Ethanol					
STORAGE TEMPERATURE	2 - 8°C					

<sup>(1)</sup> See stability table in Pre-Packed Columns Procedure for Use.

For laboratory use only. Not for use in diagnostic or therapeutic procedures.