

# TEST KITS

## Low Pressure

### CHELATING TEST KITS

DESCRIPTION	CAT. No.	COMMENTS
NICKEL CHELATE KIT Includes: 2ml HIGH Density NICKEL, 2ml LOW Density NICKEL	6BCL-KNi-2	Recommended for proteins that are easy to separate
NICKEL CHELATE KIT + 20 empty mini columns	6BCL-KNiSP-2	
NICKEL & COBALT CHELATE KIT Includes: 2ml HIGH Density NICKEL, 2ml LOW Density NICKEL, 2ml HIGH Density COBALT	6BCL-KNiCo-2	Recommended for proteins that are easy to separate
NICKEL & COBALT CHELATE KIT+ 30 empty mini columns	6BCL-KNiCoSP-2	
ZINC CHELATE KIT Includes: 2ml HIGH Density ZINC, 2ml LOW Density ZINC	6BCL-KZn-2	Recommended for proteins that are difficult to separate
ZINC CHELATE KIT + 20 empty mini columns	6BCL-KZnSP-2	
ZINC & COPPER CHELATE KIT Includes: 2ml HIGH Density ZINC, 2ml LOW Density ZINC, 2ml LOW Density COPPER	6BCL-KZnCu-2	Recommended for proteins that are difficult to separate
ZINC & COPPER CHELATE KIT + 30 empty mini columns	6BCL-KZnCuSP-2	
HIGH DENSITY CHELATE KIT Includes: 2ml HIGH Density METAL FREE, 2ml HIGH Density NICKEL, 2ml HIGH Density ZINC, 2ml HIGH Density COBALT	6BCL-KH-2	High Binding Capacity
HIGH DENSITY CHELATE KIT + 40 empty mini columns	6BCL-KHSP-2	
LOW DENSITY CHELATE KIT Includes: 2ml LOW Density METAL FREE, 2ml LOW Density NICKEL, 2ml LOW Density ZINC, 2ml LOW Density COPPER	6BCL-KL-2	High Selectivity
LOW DENSITY CHELATE KIT + 40 empty mini columns	6BCL-KLSP-2	

### PROTEIN A TEST KIT

DESCRIPTION	CAT. No.	PACK SIZE
PROTEIN A Test Kit	PA09-K-01	1 x 100 µl

### GLYOXAL TEST KIT

DESCRIPTION	CAT. No.	COMMENTS
GLYOXAL KIT Includes: 2ml LOW Density GLYOXAL 4BCL, 2ml HIGH Density GLYOXAL 4BCL, 2ml LOW Density GLYOXAL 6BCL, 2ml HIGH Density GLYOXAL 6BCL, 2ml VERY HIGH Density GLYOXAL 6BCL	GLYOXK-2	Immobilization for basic groups (Lys)

### AMINOETHYL TEST KIT

DESCRIPTION	CAT. No.	COMMENTS
AMINOETHYL LOW KIT Includes: 2ml VERY LOW Density AMINOETHYL 4BCL, 2ml LOW Density AMINOETHYL 6BCL	AMINOL-2	Immobilization for acidic groups (Asp, Glu)
AMINOETHYL COMPLETE KIT Includes: 2ml VERY LOW Density AMINOETHYL 4BCL, 2ml LOW Density AMINOETHYL 6BCL, 2ml HIGH Density AMINOETHYL 4BCL, 2ml HIGH Density AMINOETHYL 6BCL	AMINOC-2	Immobilization for acidic groups (Asp, Glu)